Abstract Submitted for the DAMOP11 Meeting of The American Physical Society

Techniques for super-resolution microscopy using NV-diamond ALEXEI TRIFONOV, Harvard, DAVID GLENN, NIR BAR-GILL, DAVID LE SAGE, RONALD WALSWORTH, Harvard-Smithsonian — We discuss the development and application of techniques for super-resolution microscopy using NV centers in diamond: stimulated emission depletion (STED), metastable ground state depletion (GSD), and stochastic optical reconstruction microscopy (STORM). NV centers do not bleach under optical excitation, are not biotoxic, and have long-lived electronic spin coherence and spin-state-dependent fluorescence. Thus NV-diamond has great potential as a fluorescent biomarker and as a magnetic biosensor.

David Le Sage Harvard-Smithsonian

Date submitted: 03 Feb 2011 Electronic form version 1.4