Abstract Submitted for the MAR08 Meeting of The American Physical Society

Direct detection of force gradient using atomic force microscopy with very small oscillation amplitude TOSHU AN¹, ATSUSHI NOMURA, TAKAHIRO NISHIO, TOYOAKI EGUCHI, KOTONE AKIYAMA², YUKIO HASEGAWA, The Institute for Solid State Physics, The University of Tokyo — Using a quartz, which has a self-sensing capability, simple configuration of atomic force microscopy (AFM) is realized, and because of its high stiffness frequencymodulation (FM) operations with a small oscillation amplitude below 100 pm is possible. The small amplitude AFM operation enhances sensitivity of short range forces. Moreover, a force gradient can be directly detected from force vs. distance measurements. We carried out FM-AFM at low temperature (LT) in ultra-high vacuum using a quartz length-extension resonator with a tungsten tip (Spring constant and resonant frequency of the resonator are 540 000 N/m and 1 MHz, respectively) (An et al., APL 76, 133114 '05). The system was simply made by attaching an extra electrode to the tip-holder of our LT-scanning tunneling microscope cooled by ³He (APL 88, 113115 '06). Direct detection of the force gradient was performed on the Si(111) 7x7 surface using very small oscillation amplitude of 70 pm at 2.4 K.

¹PRESTO, Japan Science and Technology Agency ²present address, Institute for Materials Research, Tohoku Univ.

Toshu An The Institute for Solid State Physics, The University of Tokyo

Date submitted: 12 Dec 2007

Electronic form version 1.4