

Abstract Submitted
for the MAR09 Meeting of
The American Physical Society

Hydrogen absorption by a Pd film detected by microgravimetry¹

J.I. AVILA, R. TRABOL, U.G. VOLKMANN, A.L. CABRERA, Pontificia Universidad Catolica de Chile, C. ROMERO, P. LIEVENS, Katholieke Universiteit Leuven — A thin film of palladium (200 Å) was e-beam evaporated onto a quartz crystal used in a SQM-160 Microbalance from Sigma Instruments. The crystal was then successively exposed to different hydrogen pressures, in situ, and the change in the resonant frequency was recorded as a function of time. Hydrogen absorption by the Pd film can be obtained by the shift of the resonant frequency upon hydrogenation at different pressures. A 14 Hz shift is obtained when the Pd film is saturated with hydrogen at 10 Torr of pressure. The different faces of the Pd-H system are observed with this technique.

¹Chilean grants FONDECYT 1060634, 1060628, Fund for Scientific Research, Flanders (FWO-Vlaanderen), Flemish Concerted Research Actions (GOA) and Belgian Federal Interuniversity Poles of Attraction (IAP) programs.

Alejandro Cabrera
Pontificia Universidad Catolica de Chile

Date submitted: 17 Dec 2008

Electronic form version 1.4