

Abstract Submitted
for the MAR12 Meeting of
The American Physical Society

Looking Down at Adsorption Dynamics: Filament Height Measurements with TIRFM DAVID WELCH, ZVONIMIR DOGIC, Brandeis University — Polymer adsorption is important in several contexts such as chromatography, colloid stabilization, and bio-fouling. Despite a cross-disciplinary interest in the subject, there are not many techniques to observe single-molecule adsorption events in real time. We use TIRF microscopy to accomplish this using the biological polymer f-actin adsorbed to a microscope slide via the well-known depletion interaction. We find TIRFM is able to quantitatively measure filament height, and we compare our results to theoretical predictions and previous results obtained from other systems.

David Welch
Brandeis University

Date submitted: 11 Nov 2011

Electronic form version 1.4