

Abstract Submitted
for the MAR16 Meeting of
The American Physical Society

Capturing high temperature protein conformations for low-temperature study using ultra-fast cooling DAVID MOREAU, HAKAN ATAKISI, ROBERT THORNE, Cornell University — protocols for cooling biomolecular crystals for x-ray cryocrystallography are poorly controlled, leading to crystal-to-crystal and within-crystal non-isomorphism. Furthermore, cooling times below the protein-solvent glass transition of .1 s provide ample time for biological temperature conformations to depopulate and shift. To address these issues, methods and apparatus for cooling biomolecular crystals at rates approaching 100,000 K/s have been developed. These cooling rates are sufficient to eliminate ice formation on cooling without use of cryoprotectants, and to quench additional high-temperature conformations for low-temperature study. Time scales for conformational relaxation can be characterized using variable cooling rates. Possible extension of these methods to maximize conformational quenching will be discussed.

David Moreau
Cornell University

Date submitted: 24 Nov 2015

Electronic form version 1.4