Abstract Submitted for the DPP10 Meeting of The American Physical Society

Measurement of Tearing Mode Structure and Evolution Through Fluctuation Analysis of the Motional Stark Effect Diagnostic¹ J.D. KING, M.A. MAKOWSKI, C.T. HOLCOMB, S.L. ALLEN, D.N. HILL, W.H. MEYER, R. GEER, LLNL, M.A. VAN ZEELAND, GA, T.L. RHODES, UCLA, E.C. MORSE, UC Berkeley — Fluctuation analysis of the motional Stark effect (MSE) diagnostic allows for the direct internal measurement of \tilde{B}_z . Coherent fluctuations from internal tearing modes are easily observed on several MSE channels. At the heart of this measurement is the accurate recovery of MHD sidebands at 2nd harmonic photoelastic modulator (PEM) frequencies plus/minus the MHD toroidal mode frequency of rotation. Single value decomposition (SVD) filtering and analysis has been applied for adjacent MSE channels improving signal-to-noise and enabling successful recovery of sidebands. A sample temporal and spatially resolved \tilde{B}_z measurement of a 2/1 tearing mode is presented.

¹Work performed under the auspices of the US Department of Energy under DE-AC52-07NA27344, DE-FC02-04ER54698, DE-FG03-08ER54984, and DE-FG03-89ER51116.

Josh King Lawrence Livermore National Laboratory

Date submitted: 16 Jul 2010 Electronic form version 1.4