

Abstract Submitted
for the MAS15 Meeting of
The American Physical Society

Microwave Assisted Helicon Plasmas.¹ JOHN MCKEE, UMAIR SIDDIQUI, MIGUEL HENRIQUEZ, ZACH SHORT, EARL SCIME, West Virginia University — Up to 1.2 kW of pulsed 2.45 GHz microwaves are injected into argon and helium helicon plasmas at 6 to 20 mTorr neutral pressure, at 500 W of continuous rf power, and up to 1 kG magnetic field strengths. The objective is to heat the tail of the electron energy distribution function (EEDF) and populate ion metastable states for investigation with laser-induced-fluorescence. Langmuir probes are used to measure the EEDF and optical emission spectroscopy is used to monitor ion emission from excited states populated by the additional microwave power. The injection of microwave power in argon helicon plasmas is shown to heat the high energy tail of the EEDF without increasing the plasma density. Argon ion emission is shown to increase by a factor of 4. Injection of microwaves into a helium helicon plasma is shown to cool the bulk of the EEDF and increase the plasma density. Previously absent helium ion emission lines are observed with the injection of microwaves.

¹This work is supported by US National Science Foundation grant number PHY-1360278.

John McKee
West Virginia University

Date submitted: 29 Sep 2015

Electronic form version 1.4