

Abstract Submitted
for the MAR07 Meeting of
The American Physical Society

Sorting Category: 06.8 (E)

High Frequency ESR study of the two-dimensional Heisenberg antiferromagnet copper pyrazine perchlorate. JOHAN VAN TOL, SARITHA NELLUTLA, National High Magnetic Field Lab/Florida State University, CHRIS LANDEE, Clark University, Department of Physics — Copper pyrazine perchlorate, $\text{Cu}(\text{Pyrazine})_2(\text{ClO}_4)_2$, is thought to be an almost perfect 2D Heisenberg antiferromagnet. We have performed single crystal electron magnetic resonance experiments at 120, 240, and 336 GHz on this system. Both the line-width and resonance shift are field-dependant and indicate that the 3D antiferromagnetic ordering temperature increases from 4.2 K at zero field to about 7 K at 12 Tesla. This can be interpreted by a the field-induced XY-behavior, as has been predicted[1]. A. Cuccoli et al., *Phys. Rev. B* **68**, 060402 (2003).

Prefer Oral Session
 Prefer Poster Session

Johan van Tol
vantol@magnet.fsu.edu
National High Magnetic Field Lab/Florida State University

Date submitted: 20 Nov 2006

Electronic form version 1.4